

The Münster TPC – design, DAQ and processing of a dual-phase TPC

*Institut für Kernphysik
Münster,
October 10th, 2015*

Axel Buß

Westfälische Wilhelms-Universität Münster


Email: a_buss09@uni-muenster.de



- Functionality of a dual phase TPC
- Setup of the Münster TPC
- Calibration of the photomultiplier tubes (PMTs)
- Conclusion and outlook


Motivation of direct dark matter detection and to build an institutional TPC

- WIMP (weakly interacting massive particle) as a candidate for dark matter
- Research on direct dark matter detection in various experiments (XENON, LUX, CDMS, EDELWEISS, CRESST, ...)





Functionality of a dual phase time projection chamber (TPC)

- Liquid and gaseous xenon in a cylindrical detector
→ Dual phase TPC
- Arrays of photomultipliers (PMTs) at top and bottom
- Surface of liquid xenon 1-2 centimeters below top PMTs





Functionality of a dual phase TPC

- Scintillation on electronic or nuclear recoil
→ S1 signal (short)
- More efficient light collection at the bottom due to (total) reflection at the surface of the liquid xenon
- PTFE (Teflon™) walls with high reflection properties in the VUV spectrum (178nm scintillation wavelength)





Creation of secondary signal S2

- Ionization of Xe emits electrons
- Electrons drift with constant velocity of $\sim 1.7\text{mm}/\mu\text{s}$ @ $1\text{kV}/\text{cm}$ towards surface due to high electric field
- Cathode potential: $\sim -17\text{kV}$
- Drift time \rightarrow z position




Creation of secondary signal S2

- Within gaseous phase secondary scintillation of Xe due to stronger electric field (~ 10 kV/cm) and accelerated electrons
 \rightarrow S2 signal (wide)
- Hit pattern \rightarrow xy position
- 3D position reconstruction
- Anode potential:
 $\sim +4$ kV




The Münster time projection chamber (TPC)

- Internal components
 - 14 photo multiplier tubes
 - glass fiber for LED calibration
- Dimensions:
 - inner diameter: 80 mm
 - maximum drift length: 170 mm
 - Total xenon amount: ~3kg




The Münster time projection chamber (TPC)

- Internal components
 - 14 photo multiplier tubes
 - glass fiber for LED calibration
- Dimensions:
 - inner diameter: 80 mm
 - maximum drift length: 170 mm
 - Total xenon amount: ~3kg



The Münster time projection chamber (TPC)

- Internal components
 - 14 photo multiplier tubes
 - glass fiber for LED calibration
- Dimensions:
 - inner diameter: 80 mm
 - maximum drift length: 170 mm
 - Total xenon amount: ~3kg



Photomultiplier tube (PMT)


- Emission of electrons due to the photoelectric effect
- Bialkali compounds as cathode material for high quantum efficiency (30% at 175nm)

$$\text{Q.E.} = \frac{n_{\text{detected photons}}}{n_{\text{incoming photons}}}$$


- Multiple dynodes (10) with high secondary electron emission coefficient (gain)

$$\text{gain} = \frac{N \cdot e_{\text{detected}}}{e_{\text{photo}}} \approx 2 \cdot 10^6 \text{ @ } 800 \text{ V}$$

- High time resolution ($O(1\text{ns})$)




Gain calibration





Gain calibration

Histogram of the electron output



Credit: Julian Blanke



Credit: Julian Blanke

Conclusion and outlook

- Experimental TPC in Münster for developing new techniques large TPCs might benefit from
- Gain calibration of PMTs at room temperature with external LED


Outlook:

- Refill the TPC after upgrade
- Measurements and analysis:
- Investigate S2/S1 signal ratio with different drift fields
 - 3D reconstruction
 - Characterize detector and compare it to Monte Carlo simulation (GEANT4)
 - Internal and external energy calibration

Thanks for your attention

Data Acquisition (DAQ)


- Data Acquisition by 2 CAEN1724© flash ADCs with FPPDAQ
- Trigger via external function generator or internal threshold
- Zero length encoding (ZLE) limits the amount of raw data
- Data structure:
 - Event: contains for each channel one waveform
 - Waveform: list of ADC samples




Scintillation of xenon

- Excitation or ionization of Xe atoms

Excitation:




Ionization:




Data Acquisition (DAQ)

- Data Acquisition by 2 CAEN1724© flash ADCs with FPPDAQ
- Trigger via external function generator or internal threshold
- Zero length encoding (ZLE) limits the amount of raw data
- Data structure:
 - Event: contains for each channel one waveform
 - Waveform: list of ADC samples




Hitfinding

- Hitfinding: a different approach as in Xerawdp
- Hit: signal over a given threshold of about 4-8 sigma above baseline for each individual channel





Clustering and Peak processing


- Coincidental hits are classified into peaks via a clustering plug-in
 - Peak as a sum of all contributing hits (witness)
- Clustering is followed by calculation of peak properties (>20):
 - Height
 - Area
 - Left and right boundary
 - n_contributing_channels



Distinction of electronic or nuclear recoil signal

- Excitation or Ionization dependent on primary signal
 - Electronic recoil:
 - Relatively more Ionization → High S₂/S₁ signal ratio
 - Nuclear recoil:
 - Less Ionization → Low S₂/S₁ signal ratio
- Particle identification





Properties of xenon

- Xenon as liquid scintillator (178nm wavelength)
- Heaviest non-radioactive noble gas (atomic number: 54)
 - WIMP cross section increases with A^2
 - High self-shielding
- Low activity of radioactive isotopes
- Relatively equal abundance of even and odd isotopes
 - detection of spin-dependent interactions possible
- Easy cryogenics (-100° C)
- As a liquid easily scalable

